

Märzhäuser Wetzlar GmbH & Co. KG

In der Murch 15 | 35579 Wetzlar | Germany Tel.: +49 6441 9116-0 | Fax: +49 6441 9116-40 info@marzhauser.com | www.marzhauser.com







Slide Express 2

Specifications Overall System

Automated Sample Handling for Microscopy.

	120 slides	60 slides
	76 × 26 × 1 mm	76 × 52 × 1 mm
Capacity	in 3 magazines with 20 clips each	in 3 magazines with 20 clips each
	with 2 slides each	with 1 slides each
	2 5	man i siides eddi.
Changeover time 1	approx. 15 seconds (per slide)	
Compatibility	Zeiss Axio Imager.Z2, Olympus BX61/63, Nikon Eclipse Ni-E, others on request	
Footprint	97.1 × 51.2 cm (W × L)	

Weight	handling system (without magazines) scanning stage baseplate with adapation 1 magazine incl. clips and slides ¹	approx. 65 kg approx. 6 kg approx. 26 kg approx. 2.5 kg
--------	---	--

approx. 102 cm

Specifications Scanning Stage

Accuracy ²	±1 μm
Repeatability ²	< 1 µm (bidirectional)
Resolution	0.05 μm (smallest step size)
Travel speed	max. 100 mm/s (with 4 mm ball screw pitch)

Order Information

Order information can be found overleaf.

Product Features

- integrated SCAN^{plus} scanning stage
- innovative clip concept enables high process reliability and high repeatability
- integrated measuring systems ensure position monitoring in all axes for safe sample handling
- vibration damped feet absorb vibrations and prevent slipping of the system
- quick and easy set up and commissioning
- low maintenance, high throughput
- Liquid Dispenser for application of immersion medium available as an optional accessory
- prepared for barcode scanner integration

For more accessories for microscopy please visit: www.marzhauser.com



Total height


² measured according to VDI/DGQ 3441

 $^{^{1}}$ when using clips for 2 slides 76 \times 26 mm



Märzhäuser Wetzlar GmbH & Co. KG

In der Murch 15 | 35579 Wetzlar | Germany Tel.: +49 6441 9116-0 | Fax: +49 6441 9116-40 info@marzhauser.com | www.marzhauser.com



SlideExpress 2	Order No. ¹
Handling system, scanning stage, baseplate with adapation (without magazines)	XX-24-905-0000

Magazines	Order No.
Magazine incl. 20 clips for 2 slides 76×26 mm	00-24-905-3021
Magazine incl. 20 clips for 1 slide 76 \times 52 mm	00-24-905-3022
Magazine without clips	00-24-904-3810

Clips ²	Order No. 20 Pieces	Order No. 5 Pieces
Clip for 2 slides 76 × 26 mm	00-24-905-0821	00-24-905-0851
Clip for 1 slide 76 × 52 mm	00-24-905-0822	00-24-905-0852



¹ order no. depends on microscope type ² single clips available as spare parts